



GLP-473/0.6~5mW


**BLUE LASER POINTER
AT 473nm**




GLP-473 is the first blue laser pointer in the world. It is novel, unique and amazing.



SPECIFICATIONS

Wavelength (nm)	473±1
Operating mode	CW
Output power (mW)	<5
Transverse mode	Near TEM ₀₀
Beam diameter at the aperture (1/e ² , mm)	~1.5
Beam divergence, full angle (mrad)	<1.5
Switch	Push button
Surface color	Blue
Power Consumption	DC 3V
Power supply	1×CR2 3V lithium batteries
Expected lifetime (hours)	5000
Warranty	6 months
Remark	PGL-III-473 portable laser (FDA approval) with key switch up to 50mW is available.



GLP-473	Dimensions	Accessory
 <p>Φ16.5×150(L) mm², 0.1kg</p>		 <p>1×CR2/3V lithium battery</p>